



Eeva Test: The GPS of Time-Lapse

2015 ABB CRB Workshop

Shehua Shen, MD, ELD (ABB)

Vice President, Medical & Scientific Affairs, Auxogyn

Challenges of Time Lapse

1. Time Lapse = Lots of information

- True
- False
- Don't know




2. Lots of information = useful information

- True
- False
- Don't know

3. Do you know what information Time Lapse gives to you?

- Yes
- No
- Maybe

Abnormal Cleavage




| | Blast Rate | Impl Rate |
|-----------------|------------|-----------|
| Control (n=524) | 43% | 18% |
| With AC (n=115) | 12% | 4% |
| p-value | <0.0001 | 0.05 |

- AC1 and AC2 embryos are often selected for Day 3 transfer (28.6%)
- AC embryos are often good quality (46.9% 6-10 cells, $\leq 10\%$ frag)
- Morphology is unable to detect AC embryos
- Implantation Rate: 3.7%




Athayde Wirka, et al., Atypical embryo phenotypes identified by time-lapse microscopy: high prevalence and association with embryo development, *Fertility & Sterility*. 101(6):1637-48, (2014)

Abnormal Syngamy

Normal Syngamy



Abnormal Syngamy (AS)






| | Blast Rate | Impl Rate |
|-----------------|------------|-----------|
| Control (n=443) | 45% | 18% |
| With AS (n=163) | 22% | 0% |
| p-value | <0.0001 | 0.08 |

- AS is associated with poorer developmental potential
- Many AS embryos have good morphology on Day 3 and Day 5 and are selected for transfer or freezing
- AS may be related to centrosomes from abnormal sperm

Athayde Wirka , et al., Atypical embryo phenotypes identified by time-lapse microscopy: high prevalence and association with embryo development, *Fertility & Sterility*. 101(6):1637-48, (2014)

Abnormal First Cytokinesis (A1^{cyt})



| | Blast Rate | Impl Rate |
|--------------------------------|------------|-----------|
| Control (n=443) | 45% | 17% |
| With A1 ^{cyt} (n=196) | 22% | 6% |
| p-value | <0.0001 | 0.1 |

- A1^{cyt} phenotype is associated with poorer developmental potential
- Previously research has correlated 1st cytokinesis timing (P1) to developmental competence [1]
- Combining A1^{cyt} phenotype and P1 timing may more finely discriminate embryos for de-selection

Athayde Wirka, et al., Atypical embryo phenotypes identified by time-lapse microscopy: high prevalence and association with embryo development, *Fertility & Sterility*. 101(6):1637-48, (2014). [1] Wong et al. *Nature Biotechnology* (2010)


Biological Parameters

Time-lapse observations:

- Abnormal cleavage
- Reverse cleavage
- Multinucleation
- Fragmentation dynamics
- Blastocyst collapsing and re-expansion
- ...

Time-lapse markers

Reviewed by Kaser and Racowsky, HRU 2014



More Challenges of Time Lapse

1. Do you know what “algorithm” means and what “statistical modeling” is?

- Yes
- No
- Don't know

2. How much time do you have to grade embryos and prepare for embryo transfer for each case, on an average basis?

- 15 minutes
- 30 minutes
- Unlimited time

Time-lapse parameters: Just watching is NOT enough



Two Examples:

1. We need to watch and *see*
2. We need to *see beyond* what human vision allows


Example 1: We need to watch and *see*



Example 2: We need to see beyond



Example 2: We need to see beyond




The Eeva Test

Biological
Parameters

Statistical
Modeling

Automation
+ Computer
Vision

Clinical
Validation



Regulatory
Clearance


Proven
benefit to
patients

The Eeva Test is **prognostic** for embryo selection

The Eeva Test provides **objective** information about the developmental potential of embryos



STANFORD
UNIVERSITY



nature
biotechnology

Non-invasive imaging of human embryos before embryonic genome activation predicts development to the blastocyst stage

Connie C Wong^{1,2,7}, Kevin E Loewke^{1-3,6,7}, Nancy L Bossert⁴, Barry Behr², Christopher J De Jonge⁴, Thomas M Baer⁵ & Renee A Reijo Pera^{1,2}



Top 10 Medical Breakthroughs of 2010

TOP 10 MEDICAL BREAKTHROUGHS

8. Predicting IVF Success

By Alice Park | Thursday, Dec. 09, 2010

For couples choosing to start a family with in vitro fertilization (IVF), the odds are not always in their favor. The procedure, even under the best circumstances, has a 30% chance of resulting in a live birth on average. So it was welcome news indeed when Stanford University researchers reported on a new method for selecting the strongest embryos, which would most likely result in a pregnancy and live birth.

Inside the Eeva System...


Identify Cell Divisions



Calculate Timing Intervals




Classify Embryos



P2: 10 hrs. 10 min
P3: 10 min

HIGH



P2: 14 hrs. 15 min
P3: 10 min

LOW

Day 3 Demo 1: 38 year old with 6 embryos




| Well | D3 cell number | D3 symmetry | D3 fragmentation (%) | | | Fate |
|------|----------------|-------------|----------------------|--|--|------|
| A1 | 6 | Moderate | 1-10 | | | |
| A2 | 8 | Symmetric | 1-10 | | | |
| B1 | 6 | Moderate | 1-10 | | | |
| B2 | 6 | Moderate | 1-10 | | | |
| C1 | 8 | Symmetric | 1-10 | | | |
| C2 | 8 | Symmetric | 1-10 | | | |

Desired SET

Day 3 Demo 1: 38 year old with 6 embryos



| Well | D3 cell number | D3 symmetry | D3 frag |
|------|----------------|-------------|---------|
| A1 | 6 | Moderate | |
| A2 | 8 | Symmetric | |
| B1 | 6 | Moderate | |
| B2 | 6 | Moderate | |
| C1 | 8 | Symmetric | |
| C2 | 8 | Symmetric | |



Desired SET

Day 3 Demo 1: 38 year old with 6 embryos



| Well | D3 cell number | D3 symmetry | D3 fragmentation (%) | Eeva Result | Notes | Fate |
|------|----------------|-------------|----------------------|-------------|-------|-------------|
| A1 | 6 | Moderate | 1-10 | High | | |
| A2 | 8 | Symmetric | 1-10 | Low | AC2 | |
| B1 | 6 | Moderate | 1-10 | Low | | |
| B2 | 6 | Moderate | 1-10 | Low | | |
| C1 | 8 | Symmetric | 1-10 | High | | Transferred |
| C2 | 8 | Symmetric | 1-10 | Low | AC1 | |

Outcome: Clinical pregnancy (singleton)

Day 3 Demo 4: 26 year old with 11 embryos

Fate


| Well | D3 cell number | D3 symmetry | D3 fragmentation (%) | | | Fate |
|------|----------------|-------------|----------------------|--|--|------|
| A1 | 8 | Symmetric | 1-10 | | | |
| A2 | 8 | Symmetric | 1-10 | | | |
| A3 | 8 | Symmetric | 1-10 | | | |
| B1 | Morula | Symmetric | 1-10 | | | |
| B2 | Morula | Symmetric | 1-10 | | | |
| B3 | 8 | Symmetric | 1-10 | | | |
| C1 | 8 | Symmetric | 1-10 | | | |
| C2 | 8 | Symmetric | 1-10 | | | |
| C3 | 8 | Symmetric | 1-10 | | | |
| D1 | 8 | Symmetric | 1-10 | | | |
| D2 | 6 | Symmetric | 1-10 | | | |

Desired DET

Day 3 Demo 4: 26 year old with 11 embryos



| Well | D3 cell number | D3 symmetry |
|------|----------------|-------------|
| A1 | 8 | Symmetric |
| A2 | 8 | Symmetric |
| A3 | 8 | Symmetric |
| B1 | Morula | Symmetric |
| B2 | Morula | Symmetric |
| B3 | 8 | Symmetric |
| C1 | 8 | Symmetric |
| C2 | 8 | Symmetric |
| C3 | 8 | Symmetric |
| D1 | 8 | Symmetric |
| D2 | 6 | Symmetric |



Day 3 Demo 4: 26 year old with 11 embryos



| Well | D3 cell number | D3 symmetry | D3 fragmentation (%) | Eeva Result | Notes | Fate |
|------|----------------|-------------|----------------------|-------------|-------|-------------|
| A1 | 8 | Symmetric | 1-10 | High | | Transferred |
| A2 | 8 | Symmetric | 1-10 | High | | Transferred |
| A3 | 8 | Symmetric | 1-10 | Low | | |
| B1 | Morula | Symmetric | 1-10 | Low | | |
| B2 | Morula | Symmetric | 1-10 | Low | RC | |
| B3 | 8 | Symmetric | 1-10 | High | | |
| C1 | 8 | Symmetric | 1-10 | Low | | |
| C2 | 8 | Symmetric | 1-10 | Low | RC | |
| C3 | 8 | Symmetric | 1-10 | Low | | |
| D1 | 8 | Symmetric | 1-10 | Low | | |
| D2 | 6 | Symmetric | 1-10 | Low | | |

Outcome: Clinical pregnancy (twins)

Scientific foundation of the Eeva Test

Cell division time-intervals (“P1, P2, P3”)



Recently examined for aneuploidy⁵

Correlate to implantation & blastocyst quality²⁻⁴

Reflect underlying molecular health²


Provide distinct timing windows¹

Predict successful development to blastocyst



¹Wong et al. *Nature Biotechnology* (2010), ²Meseguer et al. *Human Reprod* (2011), ³Hashimoto et al. *Fertility & Sterility* (2012), ⁴Cruz et al. *RBM Online* (2012), ⁵Chavez et al. *Nature Communications* (2012)

Summary of Early Predictive Time-Lapse Markers




* 5-8 cell stage: 40 ± 10 m; 8 cell stage: 23 ± 1 h; 9-16 cell stage: 55 ± 15 m

** dynamic assessment of fragmentation was also included in the study

Chen et al. *Fertility & Sterility*, (2013)

Ongoing Effort to Continue Scientific Discovery

- Publications and abstracts continue to increase
- Impact of Eeva Test on clinical pregnancy and implantation is under study



Publications

Clinical Validation and Research

Using the Eeva™ Test adjunctively to traditional day 3 morphology is informative for consistent embryo assessment within a panel of embryologists with diverse experience

M.P. Diamond, V. Suraj, E.J. Behnke, X. Yang, M.J. Angle, J.C. Lambe Steinmiller, R. Watterson, K. Athayde Wirka, A.A. Chen, S. Shen, *Journal of Assisted Reproduction and Genetics*, 2014

Atypical embryo phenotypes identified by time-lapse microscopy: high prevalence and association with embryo development

K. Athayde Wirka, A.A. Chen, J. Conaghan, K. Ivani, M. Gvakharia, B. Behr, V. Suraj, L. Tan, S. Shen, *Fertility and Sterility*, 2014.
[http://www.fertstert.org/article/S0015-0282\(14\)00203-9/abstract](http://www.fertstert.org/article/S0015-0282(14)00203-9/abstract)

Improving embryo selection using computer-automated time-lapse imaging plus day 3 morphology: results from a prospective multi-center trial

J. Conaghan, A.A. Chen, S.P. Willman, Kivani, P.E. Chenette, R. Boostanfar, V.L. Baker, G.D. Adamson, M.E. Abusief, M. Gvakharia, K.E. Loewke, S. Shen, *Fertility & Sterility*, 100(2): 412-9, 2013.
[http://www.fertstert.org/article/S0015-0282\(13\)00517-7/abstract](http://www.fertstert.org/article/S0015-0282(13)00517-7/abstract)

Automated time-lapse analysis in adjunctive use with morphology is highly informative in allowing diverse embryologists to select embryos with high developmental potential

M.P. Diamond, L. Tan, J. Conaghan, K. Ivani, A. Le, A.A. Chen, S. Shen, V. Suraj, European Society of Human Reproduction and Embryology (ESHRE) Annual Meeting 2014, Jun 29-Jul 2, 2014, Munich, Germany

The use of innovative, intelligent software and non-invasive embryo imaging to predict blastocyst formation by day 3

H. Marsden, J. Conaghan, K. Ivani, R. Gregoire, C. Kingsland, S. Troup
British Fertility Society/ Association of Clinical Embryologists (BFS/ACE) Annual Meeting 2013, Jan 3-5, 2013, Liverpool, UK
Awarded the "Prize Paper" for Clinical Achievement

Improved embryo selection accuracy using cell division characteristics defined by time-lapse and automated image analysis

A.A. Chen, K. Ivani, J. Conaghan, M. Gvakharia, A. Le, S. Shen
American Society for Reproductive Medicine (ASRM) Annual Meeting 2012, October 20-24, 2012, San Diego, California, *Fertility & Sterility* Vol.98(3):S17.

Prediction of embryo viability using validated cell division time intervals measured by time-lapse imaging

A.A. Chen, K.E. Loewke, S.P. Willman, P.E. Chenette, R. Boostanfar, S. Shen
American Society for Reproductive Medicine (ASRM) Annual Meeting 2012, October 20-24, 2012, San Diego, California, *Fertility & Sterility* Vol.98(3):S31.

Early cell cycle durations detected by time-lapse imaging predicts embryo developmental potential

S. Shen, A.A. Chen, S.P. Willman, P.E. Chenette, R. Boostanfar, V.L. Baker, M. Abusief, V. Suraj, K. Wirka, K. Loewke
European Society of Human Reproduction and Embryology (ESHRE) Annual Meeting 2012, July 1-4, 2012, Istanbul, Turkey, *Human Reproduction* Vol.27(2):ii22.

Research to Advance Science

Computer-automated time-lapse analysis test results correlate to clinical pregnancy and embryo implantation: a prospective, blinded, multi-center study

- Refer to [Auxogyn.com](http://www.auxogyn.com) for complete list of publications/abstracts

<http://www.auxogyn.com/clinical-innovation/reproductive-science-publications/>



How the Eeva System Works

Introducing the **Eeva™ System**



Simple & Easy to Use
Designed to fit into your lab workflow



Load embryos into **Eeva Dish**



1
2
3
A
B
C
D

CevaTM



Place dish onto Eeva Scope



auxogyn

DOE, JANE
ABC123

Confirm Alignment



36.9
5.8
cancel-confirm Menu auto-start


Ceva

| | |
|---|---|
| Scope 1 Press Here To Start <i>Ceva</i> | Scope 2 Press Here To Start <i>Ceva</i> |
| Scope 3 Press Here To Start <i>Ceva</i> | Scope 4 Press Here To Start <i>Ceva</i> |


Ceva



Analysis Begins



Day 1 & Day 2 Imaging



Scopes

Patients

Admin Settings

Imaging Complete

Eeva Result Generated

Patient

DOB

ID

Healthcare Professional


Export Results Report

Print Results Report

Export Patient Report

Print Patient Report

Export Videos



Day 3 Result Complete



Eeva Results Report

Page 2 of 2

Clinic Name :
IVF CLINIC

Healthcare Professional:
HENRY LEE MD

Patient Name: KLAIN, GIA
Patient DOB: 1976-03-01
Patient ID: 0123456789
Eeva Scope #: 03
Accession #: 00X-TST-001

Eeva Imaging Start: 2014-10-06, 08:00

Eeva Result Generated: 2014-10-08, 02:00

| Well # | Eeva Information | | | |
|--------|--|--|-------------|-------|
| | P2 Normal Range: 9hrs 20min - 11hrs 28min | P3 Normal Range: 0hrs 0min - 1hrs 44min | Eeva Result | Notes |
| A1 | 05h 15m | 00h 30m | LOW | |
| B1 | 02h 5m | 02h 35m | LOW | |
| C1 | 11h 00m | 00h 20m | HIGH | |
| D1 | 01h 30m | 05h 50m | LOW | |
| A2 | 10h 35m | 00h 15m | HIGH | |
| B2 | 03h 15m | 01h 25m | LOW | |
| C2 | 10h 15m | 00h 10m | HIGH | |
| D2 | 12h 55m | 00h 10m | LOW | |
| A3 | 00h 45m | 00h 30m | LOW | |
| B3 | 10h 25m | 01h 05m | HIGH | |
| C3 | 01h 20m | 02h 40m | LOW | |
| D3 | 00h 45m | 00h 20m | LOW | |

**Review Eeva Test Results
adjunctively to morphology**




Select with Confidence

First Clinically Validated Model

Consistent Test results within and across clinics

Generalizable Prediction Model



Test with Independent
Data Set



Multi-center (5 US clinics)



- 160 patients
- 1825 embryos



de novo clearance



FDA Clearance

Power to Predict required a Unique Path to Market

- The Eeva™ System was cleared through the FDA de novo process in June 2014
 - Pathway for innovative, low to moderate risk devices
 - More rigorous requirements
 - First device of its kind with prognostic assessment

| Eeva System | Other Time Lapse Systems |
|---|---|
| FDA De Novo Clearance | 510(k) Clearance |
| Assisted Reproduction Embryo Image Assessment System | Assisted Reproduction Accessories |

Which patients benefit from the Eeva Test....

While the Eeva Test can be used for any patient, here are some situations where its value is maximized:

The Eeva Test is:

- 1.** Best used when **embryo selection** is needed: patients who may have multiple good quality embryos on the day of embryo transfer (e.g. good responders, donor eggs, etc.)
- 2.** A tool to permit the embryologist to select with confidence when an **eSET** is planned.
- 3.** **Prognostic** not diagnostic. It has limited use in poor responders and patients with poor embryo quality.
- 4.** In all patients, the Eeva Test must be used as an **adjunct to morphology**, not as a substitute for a trained embryologist.



The future of the Eeva Test

Biological Parameters Are Only the Start

Time-lapse markers


Reviewed by Kaser and Racowsky, HRU 2014

Biological Parameters

Statistical Modeling

Automation + Computer Vision

Clinical Validation



Statistical Modeling for Multiple Parameters

Biological
Parameters


Statistical
Modeling

Automation
+ Computer
Vision

Clinical
Validation

**The future of the Eeva Test
will be...**

- Multiple biological parameters
- Multi-dimensional prediction space
- Novel surrogate image features extracted from videos



Benefits of Automation + Computer Vision

Biological
Parameters





Statistical
Modeling

Automation
+ Computer
Vision

Clinical
Validation

Automation
reduces the need
for manual
measurements

**Computer
vision**
detects
surrogate image
markers



Clinical Validation

Biological
Parameters

Statistical
Modeling

Automation
+ Computer
Vision


Clinical
Validation

Data Requirements for Developing a Predictive Model:

1. Prospectively collected
2. Multi-clinic and diverse, for generalizable and consistent results
3. Separate training & test datasets

Ongoing Development of the Eeva Test

Biological Parameters




Atypical embryo phenotypes identified by time-lapse microscopy: high prevalence and association with embryo development

Kelly Athayde Wirka, M.S.,^a Alice A. Chen, Ph.D.,^b Joe Conaghan, Ph.D.,^b Kristen Ivani, Ph.D.,^c Marina Gvakharia, M.D., Ph.D.,^d Barry Lehr, Ph.D.,^e Vaishali Suraj, M.S.,^f Lei Tan, Ph.D.,^g and Shehua Shen, M.D.^h

^a Auxogyn, Menlo Park; ^b Pacific Fertility Center, San Francisco; ^c Reproductive Science Center of the Bay Area, San Ramon; ^d Fertility Physicians of Northern California, Palo Alto Medical Foundation, San Jose; and ^e Stanford Fertility and Reproductive Medicine Center, Palo Alto, California

Wirka et al. *Fertil Steril* 2014

Statistical Modeling



Automation +
Computer Vision

Clinical Validation



An Analogy for Eeva



MAPS



GPS



NETWORK




Time -Lapse



Eeva



The Future : The Eeva Network



- Users in the network provide information in real-time
- Users receive real-time data analysis within their clinic and compared to all clinics
- Selection algorithms will be continuously updated to optimize results
- Better outcomes for members of the network

Path to Clinics and Patients



nature
biotechnology

Non-invasive imaging of human embryos before embryonic genome activation predicts development to the blastocyst stage

Connie C Wong^{1,2,7}, Kevin E Loewke^{1-3,6,7}, Nancy L Bossert⁴, Barry Behr³, Christopher J De Jonge⁵, Thomas M Baer⁵ & Rencor A Reijo Pera^{1,2}

Eeva Today (April 2015)

- 8 countries
- >45 clinics
- >4,000 patients tested
- >22,000 embryos imaged
- >35 publications



*Eva*TM

Select with Confidence



auxogyn